

Digit Prosthetic Components Straight and Angulated

Patient Information



Table of Contents

| | |
|--|---|
| 1. Intended use | 3 |
| 2. Warnings and Precautions | 3 |
| 3. Attachment | 3 |
| 4. Cleaning and hygiene instructions | 3 |
| 5. Maintenance..... | 3 |

1. Intended use

The Digit Prosthetic Straight/Angulated is intended to be used with the OPRA™ Implant System for both left and right thumb/finger amputations. The system is designed to be used for normal everyday activities in combination with a custom made cosmetic cover.

2. Warnings and Precautions

- Do not tighten the clamping screw without having the abutment inserted. It will deform the attachment device.
- The Digit Prosthetic Components are intended to be used together with the OPRA™ Implant System only and for thumb/finger amputations only.
- The Digit Prosthetic Components must be installed and serviced by a certified prosthetist.
- The Digit Prosthetic Components are intended to be used on one patient only.
- Not cleaning the device increase the risk for infection and cross contamination between patient and personnel.
- Not cleaning the device may cause the device to fall off.
- The components may alert metal detectors. Please ask you prosthetist for a certificate.

3. Attachment

Attach the prosthesis by tightening the Clamping screw using the supplied tool, while having the abutment inserted in the attachment part. Remove the prosthesis by detaching the screw.

4. Cleaning and hygiene instructions

The device should be cleaned on a regular basis, at least once per day. A tissue or cotton pins with soap, water or alcohol can be used.

5. Maintenance

- The Digit Prosthetic Straight/Angulated components should we visually inspected regularly for. Please contact your responsible prosthetists if you see any signs of damage, wear and fatigue.
- For any other questions, please contact your responsible prosthetist.

Integrum AB

Kroksläotts Fabriker 50

SE-431 37 MÖLNDAL

Sweden

Phone: +46 31-760 10 60

info@integrum.se

www.integrum.se

